


5 inch Display
5 Parameters
Touch Screen
Bluetooth

Comes include
PC Software & Android App

MD905
 Meditech Patient Monitors



Android App



Paper roll

Printer




PC Software

MD905 PC Software "COMPUTER" which is connect with patient monitor directly retransmit the monitoring data through web, so other observer in other place can observe and record the real time monitoring data and remote the patient monitor at the same time.

The Software which is running on the Tele-host computer, can retransmit the monitoring data to internet or other computers in LAN such as central monitoring system. This mode permit extending multi-observers in a simple way.

All observers can operate monitor's measurement module setup items, such as select ECG lead, scale, mode etc., so it make doctor can remote patient's monitor in his family via internet and start up NIBP measurement at any moment.



Visit Us: www.meditech.com.cn